


March 28, 2024



Closed-loop optimization of hydrogel formulations using dynamic light scattering

Acceleration Consortium Team:


Owen Melville, Ilya Yakavets, Yimu Zhao,
Nipun Gupta, Jeff Watchorn



Hydrogels are Awesome!



Cell culture




3D Printing



Drug Delivery

Hydrogels are Complicated!

Complex Formulations




Ref: <https://doi.org/10.1021/acs.biomac.9b01132>

Xu et al. (2020). *Biomacromolecules*

Complex Properties


Mechanical strength




Bio compatibility



Viscoelasticity



Shear-Thinning



Li et al. (2020) *Materials Science and Engineering*


Measuring Hydrogel Properties

Rheology



Slow
Difficult to Automate

Video Particle Tracking



Ref: <https://doi.org/10.1016/j.ymeth.2009.12.008>

Data intensive
Finicky

Mao et al. (2022). *Front. Bioeng. Biotechnol.*

Dynamic Light Scattering



Well plate compatible
Fast measurements

Krajina et al. (2017). *ACS Central Science*


A Hydrogel Self-Driving Lab!




1 Mixing hydrogel precursors

Inputs:


1. Precursor descriptors (substance)
2. Precursor concentrations (discrete)



2 Crosslinking



3 Dynamic Light Scattering characterization



4 Data analysis

Outputs:

1. Relaxation time
2. Gelation time
3. Viscosity
4. Viscoelasticity (G' and G'')

BO Features:


- Precursor representation (mixture constraints)
- Sequence optimization



5 Bayesian Optimization



The Future Needs Hydrogels!



Cell culture
3D Printing
Drug Delivery



Faster Data
Reliable Data
Self-Driving Labs




Healthy People
Happy Robots
Healthy Planet


The background is a dark, textured space filled with several glowing, translucent spheres of various sizes and colors, including orange, yellow, and red. A large, multi-faceted prism is positioned in the upper left corner, emitting a spectrum of light. The overall aesthetic is futuristic and scientific.

Thanks for your attention

SELF-DRIVING LABS ACCELERATING MATERIALS DISCOVERY



Various biological tissues show frequency-dependent viscoelastic features



Theory (Cai et al, 2021)

Incident light is scattered by embedded probe particles in a sample volume, and the scattering intensity is measured over time to derive the intensity autocorrelation function $g_2(\tau) = \langle I(t)I(t+\tau) \rangle \langle I(t) \rangle^{-2}$, where τ is the lag time between two time points and $\langle I(t) \rangle^{-2}$ is the normalization factor. This autocorrelation function relates to the mean-squared displacement $\Delta r^2(\tau)$ of the probe particles through the intermediate scattering function

$$g_1(\tau) = \exp\left(\frac{-q^2 \Delta r^2(\tau)}{6}\right), \quad (1)$$


where the scattering vector $q = 4\pi n \sin(\theta/2)/\lambda$. The intermediate scattering function $g_1(\tau)$ for ergodic systems is derived from

$$g_2(\tau) = 1 + (g_0 - 1) |g_1(\tau)|^2, \quad (2)$$

where $g_0 = g_2(0)$. From the mean-squared displacement, we then obtain the frequency-dependent complex modulus using the generalized Stokes-Einstein equation^{21,22}


$$G^*(\omega) = \frac{k_B T}{\pi a(i\omega) \mathcal{F}_u\{\langle \Delta r^2(t) \rangle\}} \quad (3)$$

where $k_B T$ denotes the thermal energy, a is the particle radius, and $\mathcal{F}_u\{\langle \Delta r^2(t) \rangle\}$ is the unilateral Fourier transform of the mean-squared displacement. This relation is extended from the purely viscous regime to continuum viscoelastic fluids and connects macroscopic stress relaxations to microscopic stress relaxations.^{23,35} Besides the complex modulus, another interesting and valuable interpretation of the mean-squared displacement is the creep compliance, which directly relates to the mean-squared displacement.³⁶ It should be noted that Equation 3 assumes that there is negligible impact from inertia of the probe, which is very sensitive to probe size and material properties of both the probe and the fluid. For a particle in a fluid, the timescale at which inertial effects are non-negligible can be estimated using the ratio of $m/\zeta = 2\rho_p a^2/9\eta_s$, where m is the mass of the particle, ζ is the drag coefficient, ρ_p is the density of the particle, a is the radius of the particle, and η_s is the viscosity of the fluid.³⁷ Particle sizes used in DLS μ R are on the order of 1 μ m, which requires a timescale of $\sim 6 \times 10^{-8}$ s for the inertia to dampen out in a material with a viscosity of 10^{-3} Pa·s. Thus, for the frequency range probed by DLS μ R (up to 10^6 Hz), inertial effects can be neglected.²²




DLSμ R Workflow Scenarios Summary. (Cai et al, 2021)


1 Embed tracer particles





2 Perform DLS




3 Determine Scattering Intensity Autocorrelation Function



4 Determine Complex Modulus G^*



Gelation Process of a Covalently Crosslinked PEG-based Gel.



Multi-day Rheology of a Living Composite of Human Breast Cancer Cells Encapsulated in Collagen I.

Relationship between Structure and Rheology of Hydrogels for Various Applications

Stojkov et al, 2021

Figure 5. Time sweep experiment shows the formation of the network structure of the hydrogel over time. The crossover point of G' and G'' indicated the gelation point.





Figure 7. Frequency sweep experiments are a powerful tool in rheology of hydrogels which show the crosslinking behavior as well as the presence of a reversible network.



A Hydrogel Self-Driving Lab!

